

iWMG18 - IB (US English)
Size: 306 x 150 mm (folded: 102 x 150 mm)
Material: 260gsm paper
Color: Black text on white paper (1C + 1C)
Jul. 05, 2021



iHome



POWERGLIDE Deep Tissue Percussion Massager




iWMG18 Quick Start Guide



For warranty visit:
www.ihome.com/support



What's Included:



iWMG18

USB-C Charging Cable & USB Adaptor



Travel Case

Attachments:

To attach, press the massage head firmly into place on the unit.
To remove, twist the massage head and pull up.



Ball
Overall Use On
Large and Small
Muscle Groups



U-Shape
Dual targeted muscle treatment
for knots and tense spots
including spine



Flat
The hammer-like
impact hits deeper
into muscles




Trigger
Pinpoint muscle treatment
including forearms, hands
and feet

This unit is equipped with a 15 minute auto-off feature for your safety.

Charging The Device

Fully charge battery for up to 7 hours before first use.


Connect the USB-C plug to the USB charging port on the back of the unit. Plug the adaptor into a 100-240V outlet.



Note: The battery working time depends on speed level and the pressure applied during use. The power level LEDs (4 green and 1 red) are solid when fully charged. As battery capacity diminishes, the green LEDs will turn off, one by one, from the top. When the battery is very low, the red LED will blink.

Warning: During charging, device will not operate as a safety precaution.

Controls + Indicators



How to Use

To turn the massager on, press and hold the power button for 2 seconds. The battery LED light indicators will illuminate to confirm the unit is on.

Press the power button once to begin the massage and continue to press (up to 5 times) to increase intensity. Once through intensity cycle, press again to stop the percussion motion and again to resume on low intensity. 5 LEDs represent 5 intensities.

Apply the massage head to the muscle you would like to massage.

When the massage is complete, turn the unit off by pressing and holding the power button for 2 seconds. The unit is off when all LED lights turn off.

NOTE: After 15 minutes of use, the unit will automatically shut off for your safety.

CAUTION: NEVER use this appliance on or near eyes or other highly sensitive areas.

Maintenance

TO CLEAN:

- Be sure to unplug the unit, if charging, and allow it to cool before cleaning. Use a soft, slightly damp cloth to wipe.
- NEVER allow water or any liquids to come into contact with the handle.
- This product has an internal, nonreplaceable lithium battery. This battery is not user-serviceable. Please dispose in accordance with local, state, province, and country regulations.

TO STORE:

- Store this unit in its box or a cool, dry place.
- DO NOT hang the unit by the power cord.
- The product should be stored at 25°C (77 °F).

IMPORTANT SAFETY INSTRUCTIONS

When using electrical products, especially when children are present, safety precautions must always be followed.

READ ALL INSTRUCTIONS PRIOR TO USE.

DANGER - TO REDUCE THE RISK OF ELECTRICAL SHOCK:

- ALWAYS UNPLUG THIS PRODUCT FROM THE ELECTRICAL OUTLET AFTER CHARGING AND BEFORE CLEANING.
- DO NOT reach for any electrical device that has fallen into water. Unplug it from the outlet immediately, first.
- DO NOT use while swimming, bathing or in the shower.
- DO NOT place the device in an area where it could fall or be pulled into a bathtub or sink.
- DO NOT drop or expose the device to water or other liquids.

WARNING - TO REDUCE THE RISK OF BURNS, FIRE, ELECTRIC SHOCK OR PERSONAL INJURY:

- This device should NEVER be left unattended when plugged in. Always unplug the charging cable when not in use and before removing and replacing any of the attachments.
- Close Adult supervision is required when this device is in use by, on or near children and disabled or persons with special needs.
- Use this device only for its intended purpose as described in these instructions. Only use the attachments included with this device.
- DO NOT operate this device if the cord or plug is damaged or frayed, has been dropped or damaged, or dropped into water, or is not working properly. Inspect periodically for any signs of damage. If damage is found, DO NOT USE. Return the device to a service center for examination and repair.
- Keep the cord away from any sources of heat that could damage the cord's insulation. Such as stoves, radiators, heat registers, or other appliances that produce heat.
- DO NOT operate where aerosol products are in use or where oxygen is being administered.
- DO NOT block or inhibit the proper ventilation of the device by using it under a blanket, pillow or other soft surfaces (bed or sofa). Overheating can occur and can cause personal injury, including fire or electric shock.
- ALWAYS keep the ventilation openings free of hair, lint, dust and other debris.
- ALWAYS keep long hair away from the device while in use.
- DO NOT use near loose clothing or jewelry.
- ALWAYS turn the device "Off", prior to charging the unit.
- For Indoor Use Only. DO NOT use outdoors.
- This device is designed for personal, household, nonprofessional use only.
- Only recharge the device with the supplied charger. Use of other chargers can create a risk of fire.
- Never drop or insert any object into any opening.
- Unplug unit and clean only with dry cloth. Stronger cleaning agents, such as Benzene, thinner or similar materials can damage the surface of the unit.
- The main plug of power adaptor is used as the disconnect device, it shall remain readily operable.
- Stop using immediately and seek medical attention if you experience any abnormal or severe pain.

CAUTION: IF YOU ARE PREGNANT OR HAVE AN ILLNESS, PLEASE CONSULT WITH YOUR DOCTOR BEFORE USING THIS DEVICE.

PLEASE READ ALL THE INSTRUCTIONS CAREFULLY PRIOR TO USE.

Consult with your doctor before use, if:

- You are pregnant.
- You have a pacemaker.
- You have an illness or otherwise have concerns regarding your health.
- NOT recommended for use by people with diabetes.
- NEVER leave the device unattended and plugged in around children.
- NEVER use this device directly on swollen, inflamed or open skin wounds, joints of elbows and knees, breasts or abdomen.
- THIS DEVICE IS NOT A SUBSTITUTE FOR MEDICAL ATTENTION.

INTERNAL RECHARGEABLE LITHIUM-ION BATTERY

This device uses an internal non-replaceable rechargeable lithium-ion battery that should last the lifetime of the device. If you suspect your battery may be dead, try charging it for an extended period of time and try resetting the unit.

WARNING/DISPOSAL OF UNIT

Do NOT attempt to change the battery. Doing so will void the warranty and could cause injury and damage to your device. To dispose of the unit, follow the instructions below.

- Never dispose of used batteries/electronics with ordinary solid wastes, since they may contain toxic substances.
- Always dispose of used batteries/electronics in accordance with the prevailing community regulations.
- Many electronics retailers now recycle batteries and other electronic components at no charge. To find a drop-off location near you, visit www.call2recycle.com

GENERAL BATTERY PRECAUTIONS

- Do not dispose of batteries in the household trash or fire.
- Never use or charge a battery if it appears to be leaking, discolored, deformed, or in any way abnormal.
- Never disassemble the unit or attempt to remove battery as that could cause leakage of alkaline solution or other electrolytic substances.
- Avoid charging a rechargeable battery when the temperature is not within 5°C - 35° C (41° - 95° Fahrenheit).
- Never expose a battery to any liquid.
- Never expose a battery to excessive heat such as direct sunlight or fire. Never use a battery pack in an unventilated vehicle where excessive internal temperatures may be encountered.
- Always keep batteries out of the reach of infants and small children.
- If the unit is to be unused for an extended period of time, such as a month or longer, charge the unit fully before storage.

FCC Information

This equipment complies with FCC RF radiation exposure limits set forth for an uncontrolled environment.

This device complies with Part 15 of the FCC Rules, operation is subject to the following two conditions:

- (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Warning: Changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules.

These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications.

However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

ISED Information

This equipment complies with IC RSS-102 radiation exposure limits set forth for an uncontrolled environment.

This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s). Operation is subject to the following two conditions:

- (1) This device may not cause interference.
- (2) This device must accept any interference, including interference that may cause undesired operation of the device.

The device complies with ISED RF exposure requirements, users can obtain Canadian information on RF exposure and compliance.

Any Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Cet équipement est conforme aux limites d'exposition aux radiations IC CNR-102 établies pour un environnement non contrôlé.

L'émetteur/récepteur exempt de licence contenu dans le présent appareil est conforme aux CNR d'Innovation, Sciences et Développement économique Canada applicables aux appareils radio exempts de licence. L'exploitation est autorisée aux deux conditions suivantes:

- (1) L'appareil ne doit pas produire de brouillage;
- (2) L'appareil doit accepter tout brouillage radioélectrique subi, même si le brouillage est susceptible d'en compromettre le fonctionnement.

L'appareil est conforme aux exigences d'exposition RF (SED), les utilisateurs peuvent obtenir des informations canadiennes sur l'exposition aux RF et la conformité.

Tout changement ou modification non expressément approuvé par la partie responsable de la conformité peut annuler l'autorisation de l'utilisateur à utiliser l'équipement.

CAN ICES-3(B) / NMB-3(B)
© 2021 SDI Technologies, Inc. All rights reserved.
Questions? Visit www.ihome.com

iWMG18-063021-A

Printed in China